

Rat Osteopontin ELISA Kit

Cat. No. GTX00628

实验应用 ELISA

种属反应 Rat

Package

1 kit

产品说明

Content

	Kit components	Package	Storage Temp
1.	Anti-Rat PONT ELISA micro-plate	8 wells * 12 strips	4-8°C
2.	1x Diluent Running Buffer	1 * 60 ml	4-8°C
3.	20x Concentrated Wash Solution	1 * 50 ml	4-8°C
4.	100x Detection Antibody	1 * 150 µl	4-8°C
5.	100x HRP-Sterptavidin	1 * 150 µl	4-8°C
6.	Chromogen-Substrate Solution	1 * 12 ml	4-8°C
7.	Stop Solution	1 * 12 ml	4-8°C
8.	Rat PONT Calibrator	1 vial	4 or -20°C

摘要

The osteopontin (PONT) test kit is a highly sensitive two site enzyme linked immunoassay (ELISA) for measuring PONT in biological fluid of rats.

实验应用

应用说明

Dynamic Detection Range :

4 ng/ml-0.06 ng/ml

Sensitivity :

0.02 ng/ml

属性


注意事项

仅供实验室使用。不适用于人类或动物的任何临床, 治疗或诊断用途。不适合动物或人类食用。



For full product information, images and publications, please visit our [website](#).

產品圖片



GTX00628 ELISA Image

ELISA analysis of immunogen using GTX00628 Rat Osteopontin ELISA Kit.



For full product information, images and publications, please visit our [website](#).