

## TNF alpha antibody [MT11B10] (Biotin)

Cat. No. GTX03005-02

宿主	Rat
克隆	Monoclonal
同种型	IgG2a
实验应用	ELISA, Sandwich ELISA
种属反应	Mouse

## 实验应用

## 应用说明

\*最佳稀释倍数与浓度应由研究人员确认

Suggested dilution	Recommended dilution
ELISA	Assay dependent
Sandwich ELISA	Assay dependent

**Note : Capture: GTX03007, Detection: GTX03005-02.**

以下为常规应用缩写的中文注解

WB: 免疫印迹

ICC/IF: 细胞染色

IHC-P: 石蜡切片

IHC-Fr: 冰冻切片


## 属性

形式	Liquid
存储溶液	PBS
保存剂	0.02% Sodium azide
存放说明	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
浓度	0.5 mg/ml (Please refer to the vial label for the specific concentration.)
偶联	Biotin Biotinylated through reaction with a N-hydroxysuccinimide ester of biotin.
注意事项	仅供实验室使用。不适用于人类或动物的任何临床、治疗或诊断用途。不适合动物或人类食用。



For full product information, images and publications, please visit our [website](#).

## 產品圖片

**GTX03005-02 ELISA Image**

Sandwich ELISA analysis of mouse TNF-alpha protein using GTX03007 TNF alpha antibody [MT1C8 + 23C9] as coating antibody and GTX03005-02 TNF alpha antibody [MT11B10] (Biotin) as detecting antibody.



For full product information, images and publications, please visit our [website](#).