

Interferon alpha antibody [MT1 + MT3 + MT5]

Cat. No. GTX03015

宿主	Mouse
克隆	Multidonal
同种型	IgG1
实验应用	ELISA, Sandwich ELISA
种属反应	Human, Monkey

实验应用

应用说明

*最佳稀释倍数与浓度应由研究人员确认

Suggested dilution	Recommended dilution
ELISA	Assay dependent
Sandwich ELISA	Assay dependent

Note : Capture: GTX03015, Detection: GTX03017-02.

以下为常规应用缩写的中文注解

WB: 免疫印迹

ICC/IF: 细胞染色

IHC-P: 石蜡切片

IHC-Fr: 冰冻切片


属性

形式	Liquid
存储溶液	PBS
保存剂	0.02% Sodium azide
存放说明	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
浓度	1 mg/ml (Please refer to the vial label for the specific concentration.)
偶联	Unconjugated
注意事项	仅供实验室使用。不适用于人类或动物的任何临床, 治疗或诊断用途。不适合动物或人类食用。



For full product information, images and publications, please visit our [website](#).

產品圖片



GTX03015 ELISA Image

Sandwich ELISA analysis of human IFN- α protein using GTX03015 Interferon alpha antibody [MT1 + MT3 + MT5] as coating antibody and GTX03017-02 Interferon alpha antibody [MT2 + MT4 + MT6] (Biotin) as detecting antibody.



For full product information, images and publications, please visit our [website](#).