

IL13 ELISA pair [IL13-I/IL13-3]

Cat. No. GTX03077

实验应用 ELISA, Sandwich ELISA

种属反应 Human, Monkey

Package

1 pair

产品说明

Product Content

Cat No.	Type	宿主	偶联	Package
GTX02959-02	Detection antibody	Rat	Biotin	250 µg
GTX02960	Capture antibody	Mouse	Unconjugated	250 µg

属性

形式 Liquid


存放说明

Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

注意事项

仅供实验室使用。不适用于人类或动物的任何临床、治疗或诊断用途。不适合动物或人类食用。

产品图片



GTX03077 ELISA Image

Sandwich ELISA analysis of monkey IL-13 protein using GTX02960 IL13 antibody [IL13-I] as coating antibody and GTX02959-02 IL13 antibody [IL13-3] (Biotin) as detecting antibody.



For full product information, images and publications, please visit our [website](#).

Date 2025 / 05 / 13 Page 1 of 1