

ImmunoHisto-Sealer

Cat. No. GTX03606

实验应用 IHC-P

Package
30 ml

产品说明

摘要

ImmunoHisto-Sealer is a preservation solution that acts as a sealant to improve tissue morphology and stain over an extended period of time. When applied to the stained section, specimens can be permanently mounted without chromogens fading. Additional benefits of the ImmunoHisto-Sealer include: saves hands on time when mounting slides, avoids the use of toxic chemicals such as xylene and eliminates accumulation of hazardous waste. This product is designed for tissue specimens stained with peroxidase, alkaline phosphatase-based systems and various chromogens. After treating with GTX03606 ImmunoHisto-Sealer, water solution chromogens can be mounted in Permanent Mounting Media.

属性


形式 Liquid

存放说明 Store at room temperature.

偶联 Unconjugated

注意事项 仅供实验室使用。不适用于人类或动物的任何临床，治疗或诊断用途。不适合动物或人类食用。

产品图片



GTX03606 IHC-P Image

IHC-P analysis of human tonsil tissue using Cytokeratin H.M.W antibody HRP Blue with/without GTX03606 ImmunoHisto-Sealer.



For full product information, images and publications, please visit our [website](#).