

Mouse Anti-Human IgG2 (Fd) antibody [HP6014] (Biotin)

Cat. No. GTX04181-02

宿主	Mouse
克隆	Monoclonal
同种型	IgG1
实验应用	WB, IHC-P, IHC-Fr, FCM, IP, ELISA, Depletion, Purification
种属反应	Human

Package

500 µg

实验应用

应用说明

*最佳稀释倍数与浓度应由研究人员确认

Suggested dilution	Recommended dilution
WB	Assay dependent
IHC-P	Assay dependent
IHC-Fr	Assay dependent
FCM	Assay dependent
IP	Assay dependent
ELISA	1:5000-1:10000
Depletion	Assay dependent
Purification	Assay dependent

以下为常规应用缩写的中文注解

WB: 免疫印迹

ICC/IF: 细胞染色

IHC-P: 石蜡切片

IHC-Fr: 冰冻切片

属性

形式	Liquid
存储溶液	PBS
保存剂	0.09% sodium azide
存放说明	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C.
浓度	0.5 mg/ml (Please refer to the vial label for the specific concentration.)
免疫原	Human IgG2 myeloma protein
纯化方式	Purified IgG1


For full product information, images and publications, please visit our [website](#).

Date 2025 / 07 / 03 Page 1 of 2

偶联 Biotin

注意事项 仅供实验室使用。不适用于人类或动物的任何临床，治疗或诊断用途。不适合动物或人类食用。

產品图片



GTx04181-02 ELISA Image

ELISA analysis of purified human immunoglobulins using serially diluted GTx04181-02 Mouse Anti-Human IgG2 (Fd) antibody [HP6014] (Biotin) followed by Streptavidin (HRP).



For full product information, images and publications, please visit our [website](#).