

Human TIAF1 protein, GST tag

Cat. No. GTX102565-pro**Species** Human**Package**

100 µg, 25 µg

属性

形式 Liquid**存储溶液** PBS, 0.1% SDS**保存剂** No preservatives**存放说明** Store as concentrated solution. Centrifuge briefly prior to opening vial. Aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.**浓度** Batch dependent (Please refer to the vial label for the specific concentration.)**Region/Sequence** N-terminal GST tagged full length protein (Met removed) of NP_004731**表达系统** E. coli**纯度** >90%**偶联** Unconjugated**注意事项** For In vitro laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.For full product information, images and publications, please visit our [website](#).