

## Rat Anti-Mouse IgG2a (Heavy chain) antibody [LO-MG2a-9]

Cat. No. GTX11584

|      |            |
|------|------------|
| 宿主   | Rat        |
| 克隆   | Monoclonal |
| 同种型  | IgG1       |
| 实验应用 | WB, ELISA  |
| 种属反应 | Mouse      |

Package  
250 µg

## 实验应用

## 应用说明

\*最佳稀释倍数与浓度应由研究人员确认

| Suggested dilution | Recommended dilution         |
|--------------------|------------------------------|
| WB                 | Assay dependent              |
| ELISA              | 5µg/ml (as coating antibody) |

以下为常规应用缩写的中文注解

WB: 免疫印迹

ICC/IF: 细胞染色

IHC-P: 石蜡切片

IHC-Fr: 冰冻切片

## 产品说明

This antibody recognizes the mouse gamma 2a heavy chain, binding to all allotypes and does not cross-react with other murine immunoglobulin classes or subclasses.

## 属性

|      |                                                                                                                                                                                                                            |
|------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 形式   | Liquid                                                                                                                                                                                                                     |
| 存储溶液 | PBS                                                                                                                                                                                                                        |
| 保存剂  | 0.1% Sodium azide                                                                                                                                                                                                          |
| 存放说明 | Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles. |
| 浓度   | 1 mg/ml (Please refer to the vial label for the specific concentration.)                                                                                                                                                   |
| 免疫原  | Purified IgG2a from BALB/c mice                                                                                                                                                                                            |
| 纯化方式 | Purified by affinity chromatography from tissue culture supernatant                                                                                                                                                        |
| 偶联   | Unconjugated                                                                                                                                                                                                               |
| 注意事项 | 仅供实验室使用。不适用于人类或动物的任何临床, 治疗或诊断用途。不适合动物或人类食用。                                                                                                                                                                                |

For full product information, images and publications, please visit our [website](#).