

Osteocalcin antibody [OCG4]

Cat. No. GTX13421

宿主	Mouse
克隆	Monoclonal
同种型	IgG1
实验应用	WB, IHC-P, IHC-Fr, ELISA, Sandwich ELISA
种属反应	Human, Rabbit, Goat, Sheep, Bovine, Chicken, Pig

实验应用

应用说明

*最佳稀释倍数与浓度应由研究人员确认

Suggested dilution	Recommended dilution
WB	1-10 µg/ml
IHC-P	1-10µg/ml
IHC-Fr	1-10µg/ml
ELISA	Assay dependent
Sandwich ELISA	1-10µg/ml (with immobilized antigen)

Note : Antigen retrieval with 0.4ug/ml Proteinase K.**Antigen retrieval with 0.4ug/ml Proteinase K.****Capture : GTX13421 ; Detection : GTX13419**

以下为常规应用缩写的中文注解

WB: 免疫印迹

ICC/IF: 细胞染色

IHC-P: 石蜡切片

IHC-Fr: 冰冻切片

属性

形式	Liquid
存储溶液	PBS, 1% BSA
保存剂	No preservative
存放说明	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
浓度	Batch dependent (Please refer to the vial label for the specific concentration.)
偶联	Unconjugated
注意事项	仅供实验室使用。不适用于人类或动物的任何临床、治疗或诊断用途。不适合动物或人类食用。

For full product information, images and publications, please visit our [website](#).

產品圖片

GTX13421 IHC-P Image

IHC-P analysis of various cancer tissue samples using GTX13421 Osteocalcin antibody [OCG4].



For full product information, images and publications, please visit our [website](#).