

Trans-Hydroxyproline antibody

Cat. No. GTX13469

宿主	Rabbit
克隆	Polyclonal
同种型	IgG
实验应用	ELISA
种属反应	Species independent

产品说明

摘要	Compounds	Affinity compared to Trans-hydroxyproline-BSA
	Trans-hydroxyproline-BSA	1
	Cis-hydroxyproline-BSA	1
	Proline-BSA	40 Times
	Pyrroline-5-carboxylic acid-BSA	900 Times

实验应用

应用说明

Optimal dilutions should be determined by each laboratory for each application.

属性

形式	Liquid
存储溶液	PBS, 20% Glycerol
保存剂	No preservative
存放说明	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
浓度	Batch dependent (Please refer to the vial label for the specific concentration.)
偶联	Unconjugated
注意事项	仅供实验室使用。不适用于人类或动物的任何临床，治疗或诊断用途。不适合动物或人类食用。



For full product information, images and publications, please visit our [website](#).