

## Human Ferritin Light Chain protein, His tag

Cat. No. GTX138607-pro

Species Human

Package  
100 µg

## 实验应用

Observed MW 26 kDa.

## 属性

形式 Lyophilized powder

存储溶液 Reconstitute with distilled water to 0.3 mg/ml. Lyophilized from PBS, 2% Trehalose, 5% Mannitol, 0.01% Tween-80.

保存剂 No Preservative

存放说明 Store at 4°C or below. After reconstitution, keep as concentrated solution. Aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles. For long-term storage after reconstitution, aliquot and store at -70°C or below. Do not vortex.

Region/Sequence C-terminal His tagged full-length Human Ferritin Light Chain protein (1-175 a.a. of NP\_000137.2)


表达系统 HEK293 cells

纯度 &gt;85%

偶联 Unconjugated

注意事项 仅供实验室使用。不适用于人类或动物的任何临床、治疗或诊断用途。不适合动物或人类食用。

## 產品图片



## GTX138607-pro Image

GTX138607-pro Human Ferritin Light Chain protein, His tag was analyzed using SDS-PAGE and stained with coomassie blue and captured by monochrome camera.



For full product information, images and publications, please visit our [website](#).