

USP25 blocking peptide

Cat. No. GTX17216-PEP

实验应用 Neutralizing /Inhibition

Package

50 µg

产品说明

摘要

USP25 peptide is used for blocking the activity of USP25 antibody (catalog number GTX17216). It usually blocks the antibody activity completely in Western blot by incubating the peptide with equal volume of antibody for 30 min at 37°C.

实验应用

应用说明

*最佳稀释倍数与浓度应由研究人员确认

Suggested dilution

Neutralizing /Inhibition

Recommended dilution

Assay dependent

以下为常规应用缩写的中文注解

WB: 免疫印迹

ICC/IF: 细胞染色

IHC-P: 石蜡切片

IHC-Fr: 冰冻切片

属性

形式

Liquid

存储溶液

PBS, 0.1% BSA

保存剂

0.02% Sodium azide

存放说明

Store as concentrated solution. Aliquot and store at -20°C or below. Avoid freeze-thaw cycles.

浓度

200 µg/ml (Please refer to the vial label for the specific concentration.)

偶联

Unconjugated

注意事项

仅供实验室使用。不适用于人类或动物的任何临床，治疗或诊断用途。不适合动物或人类食用。



For full product information, images and publications, please visit our [website](#).