

Mouse VCAM1 / CD106 protein, His tag

Cat. No. GTX17321-pro

Species Mouse**Package**

50 µg

属性

形式 Liquid**存储溶液** PBS, 10% Glycerol**保存剂** No preservative**存放说明** Store as concentrated solution. Aliquot and store at -20°C or below. Avoid freeze-thaw cycles.**浓度** 0.5 mg/ml (Please refer to the vial label for the specific concentration.)**Region/Sequence**

25-698aa, C-terminal His-Tag; FKIEISPEYK TIAQIGDSMA LTCSTTGCEs PLFSWRTQID SPLNAKVRTE GSKSVLTMEP
VSFENEHSYL CTATCGSGKL ERSIHVDIYS FPKDPEIQFS GPLEVGKPVV VKCLAPDIYP VYRLEIDLFK GDQLMNRQEF SSEEMTKSLE
TKSLEVTFTP VIEDIGKALV CRAKLHIDQI DSTLKERETV KELQVYISPR NTTISVHPST RLQEGGAVTM TCSSEGLPAP EIFWGRKLDN
EVLQLLSGNA TLTIAMRME DSGVYVCEGV NLIGRDKAEV ELVVQEKPI VDISPGSQVA AQVGDSVLT CAAIGCDSPS
FSWRTQTDSP LNGVVRNEGA KSTLVSSVG FEDEHSYLCA VTCLQRTLEK RTQVEVYSFP EDPVIKMSGP LVHGRPVTN
CTVPNVYFPD HLEIELLKGE TTLMKKYFLE EMGIKSLETK ILETFTIPTI EDTGKSLVCL ARLHSGEMES EPKQRQSVQP LYVNVAPKET
TIWVSPSPIL EEGSPVNLTC SSDGIPAPKI LWSRQLNNGE LQPLSENTTL TFMSTKRDDS GIYVCEGINE AGISRKSVEL IIQVSPKDIQ
LTVFPSKSVK EGDTVISCT CGNVPETWII LKKKAKTGDM VLKSVDSYIT IRQAQLQDAG IYECESKTEV GSQRLSLTLD
VKGKEHNKNY FSPELEHHHH HH

表达系统 Insect cells**偶联** Unconjugated**注意事项** For In vitro laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

For full product information, images and publications, please visit our [website](#).