

Human HMGB1 protein, His tag

Cat. No. GTX17619-pro

Species Human

Package

20 µg

产品说明

摘要

HMG1 Human Recombinant fused to an 8 aa His-Tag at C-terminus produced in High Five insect cells is a single, glycosylated, polypeptide chain (amino acids 1-215) containing 223 amino acids and having a molecular mass of 25 kDa.

属性

形式 Liquid

存储溶液 20mM Tris-HCl, 10% Glycerol, 1mM EDTA, 0.5mM DTT

保存剂 No preservative

存放说明 Store as concentrated solution. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, we recommend adding a carrier protein (0.1% HSA or BSA), then aliquot and store at -20°C or below. Avoid freeze-thaw cycles.

浓度 1 mg/ml (Please refer to the vial label for the specific concentration.)

Region/Sequence

C-terminal His-Tag; MGKGDPPKPR GKMSSYAFFV QTCREEHKKK HPDASVNFSE FSKKCSERWK TMSAKEKGKF EDMAKADKAR
YEREMKTYIP PKGETKKKFK DPNAKRPSP AFFLCSEYR PKIKGEHPGL SIGDVAKKLG EMWNNTAADD KQPYEKKA
LKEKYEKDIA AYRAKGKPD AAKGVVKA EKSKKKKEEED EDEEDEEEEE EDEEDEEEEE DDDDELEHHH HHH

表达系统 Baculovirus (High-5 cells)

纯度 > 90% by SDS-PAGE

偶联 Unconjugated

注意事项 For In vitro laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.



For full product information, images and publications, please visit our [website](#).