

## HEF1 antibody [14A11]

**Cat. No. GTX18059**

|      |            |
|------|------------|
| 宿主   | Mouse      |
| 克隆   | Monoclonal |
| 同种型  | IgG2b      |
| 实验应用 | ICC/IF, IP |
| 种属反应 | Human      |

## 实验应用

## 应用说明

**Recommended Starting Dilutions:**

at 1 $\mu$ g of antibody per 200 $\mu$ g of whole cell lysate for IP, & 1:100 for IF  
Optimal dilutions should be determined experimentally by the researcher.

## 属性

|      |  |
|------|--|
| 形式   | Liquid   |
| 存储溶液 | PBS  |
| 保存剂  | 0.1% Sodium azide  |
| 存放说明 | Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles. |
| 浓度   | Batch dependent (Please refer to the vial label for the specific concentration.)   |
| 偶联   | Unconjugated   |
| 注意事项 | 仅供实验室使用。不适用于人类或动物的任何临床, 治疗或诊断用途。不适合动物或人类食用。  |



For full product information, images and publications, please visit our [website](#).