

EPO antibody [85]

Cat. No. GTX19485

宿主	Mouse
克隆	Monoclonal
同种型	IgG1
实验应用	WB, ELISA
种属反应	Human

实验应用

应用说明

ELISA: Use at an assay dependent dilution. Optimal dilutions/concentrations should be determined by the end user.

属性

形式	Liquid
存储溶液	PBS
保存剂	No preservatives
存放说明	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
浓度	1 mg/ml (Please refer to the vial label for the specific concentration.)
偶联	Unconjugated
注意事项	仅供实验室使用。不适用于人类或动物的任何临床，治疗或诊断用途。不适合动物或人类食用。

产品图片



GTX19485 WB Image

Western blot analysis detecting EPO antigen (100 ng/well) using mouse anti-human EPO antibody (clone 85) at 1:1000 dilution (1 ug/mL)



For full product information, images and publications, please visit our [website](#).