

DcR2 / CD264 antibody

Cat. No. GTX19914

宿主	Rabbit
克隆	Polyclonal
同种型	IgG
实验应用	WB
种属反应	Human

实验应用

应用说明

WB: Use at a concentration of 2.5 - 5 µg/ml for 2 hrs at room temperature. Predicted molecular weight: 36 kDa. Optimal dilutions/concentrations should be determined by the end user.

属性

形式	Liquid
存储溶液	PBS, BSA
保存剂	Sodium azide
存放说明	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. DO NOT FREEZE.
浓度	0.5 mg/ml (Please refer to the vial label for the specific concentration.)
偶联	Unconjugated
注意事项	仅供实验室使用。不适用于人类或动物的任何临床, 治疗或诊断用途。不适合动物或人类食用。



For full product information, images and publications, please visit our [website](#).