

Factor XII Light chain antibody [C6B7]

Cat. No. GTX21008

宿主	Mouse
克隆	Monoclonal
同种型	IgG1
实验应用	Coagulation assay, WB, ELISA
种属反应	Human, Mouse

引用文献 (2)

实验应用

应用说明

*最佳稀释倍数与浓度应由研究人员确认

Suggested dilution	Recommended dilution
WB	Assay dependent
ELISA	Assay dependent
Coagulation assay	Assay dependent

以下为常规应用缩写的中文注解

WB: 免疫印迹

ICC/IF: 细胞染色

IHC-P: 石蜡切片


IHC-Fr: 冰冻切片

属性

形式	Liquid
存储溶液	PBS
保存剂	No preservative
存放说明	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
浓度	Batch dependent (Please refer to the vial label for the specific concentration.)
偶联	Unconjugated
RRID	AB_374291
注意事项	仅供实验室使用。不适用于人类或动物的任何临床, 治疗或诊断用途。不适合动物或人类食用。

For full product information, images and publications, please visit our [website](#).

產品圖片



GTx21008 WB Image

The data was published in the journal Oncotarget in 2016. [PMID: 27248165](https://pubmed.ncbi.nlm.nih.gov/27248165/)



For full product information, images and publications, please visit our [website](https://www.genetex.com).