

CD99 antibody [O13]

Cat. No. GTX23382

宿主	Mouse
克隆	Monoclonal
同种型	IgG1
实验应用	IHC-P
种属反应	Human

实验应用

应用说明

IHC-P: Use at a dilution of 1/40 or greater for Biotin-based detection systems. For optimal staining, the primary antibody should be incubated 20-60 minutes at room temperature. Sections should be pretreated with an antigen unmasking solution. Optimal dilutions/concentrations should be determined by the end user.

属性

形式	Liquid
存储溶液	PBS, 1% BSA
保存剂	0.1% Sodium azide
存放说明	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C.
浓度	Batch dependent (Please refer to the vial label for the specific concentration.)
偶联	Unconjugated
注意事项	仅供实验室使用。不适用于人类或动物的任何临床、治疗或诊断用途。不适合动物或人类食用。



For full product information, images and publications, please visit our [website](#).