

Interferon gamma antibody [H3-1]

Cat. No. GTX28318

宿主	Mouse
克隆	Monoclonal
同种型	IgG1
实验应用	WB, ELISA
种属反应	Human

实验应用

应用说明

This antibody can be used to neutralize the anti-viral activity of interferon gamma. Can also be used for the isolation and purification of recombinant interferon gamma.

属性

形式	Liquid
存储溶液	PBS
保存剂	0.1% Sodium azide
存放说明	Keep as concentrated solution, aliquot and store at 4°C. DO NOT FREEZE. We recommend to parafilm vials during storage to prevent evaporation of buffer.
浓度	Batch dependent (Please refer to the vial label for the specific concentration.)
偶联	Unconjugated
注意事项	仅供实验室使用。不适用于人类或动物的任何临床、治疗或诊断用途。不适合动物或人类食用。



For full product information, images and publications, please visit our [website](#).