

Human CXCL1 / GRO Alpha protein

Cat. No. GTX29806-pro

Species Human**Package**

25 µg

产品说明

摘要

GRO Alpha Human Recombinant produced in E.Coli is a single, non-glycosylated, polypeptide chain containing 73 amino acids and having a molecular mass of 7.811 kDa.

属性

形式 Lyophilized powder**存储溶液** Reconstitute with ddH₂O to 0.1-1.0mg/ml. Lyophilized from 20mM Phosphate buffer (pH7.4), 50mM NaCl.**保存剂** No preservative**存放说明** Store at -20°C or below. After reconstitution, keep as concentrated solution. Aliquot and avoid freeze-thaw cycles.**Region/Sequence** The sequence of the first five N-terminal amino acids was determined and was found to be Ala-Ser-Val-Ala-Thr.**表达系统** E. coli**纯度** > 98% by SDS-PAGE and HPLC.**内毒素** < 1 EU/µg**偶联** Unconjugated**注意事项** For In vitro laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

For full product information, images and publications, please visit our [website](#).