

## EpCAM antibody [ESA 214] (FITC)

Cat. No. GTX30708

宿主	Mouse
克隆	Monoclonal
同种型	IgG1
实验应用	ICC/IF, FCM, IHC
种属反应	Human

引用文献 ( 8 )

## 实验应用

## 应用说明

IHC, ICC (frozen or formalin-fixed paraffin-embedded (FFPE) tissue sections, cell smears. For IHC dilute conc. antibodies 1:50-1:100, use streptavidin~biotin system or polymer system, incubate 30 minutes at room temperature. IF, FITC~conjugated antibodies 10-20 µg/ml (1:10-1:20), incubate for 2 hours in the dark at RT or it can also be incubated overnight at 4°C. FITC~conjugated antibodies, 0.2-1.0 µg/0.1 ml (1:200-1:1,000)

## 属性

形式	Liquid
存储溶液	PBS, 1% BSA
保存剂	0.05% Sodium azide
存放说明	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. DO NOT FREEZE. Protect from light.
浓度	0.20 mg/ml (Please refer to the vial label for the specific concentration.)
偶联	Fluorescein isothiocyanate (FITC) <a href="#">Wavelength</a>
注意事项	仅供实验室使用。不适用于人类或动物的任何临床、治疗或诊断用途。不适合动物或人类食用。



For full product information, images and publications, please visit our [website](#).

Date 2025 / 07 / 21 Page 1 of 1