

## Hydroxyprogesterone antibody [HPRO-3]

Cat. No. GTX40803

宿主	Mouse
克隆	Monoclonal
同种型	IgG2b
实验应用	ELISA
种属反应	Species independent

## 产品说明

Compounds	Cross reactivity (%)
Progesterone	11
17-Hydroxyprogesterone	100
11-Hydroxyprogesterone	N/A
5-alpha-pregnane-3	N/A
Corticosterone	N/A
Pregnenolone	N/A
Deoxycorticosterone	15
Deoxycortisol	N/A
Cortisol	N/A

## 实验应用

## 应用说明

\*最佳稀释倍数与浓度应由研究人员确认

Suggested dilution	Recommended dilution
ELISA	Assay dependent

**Note : used in competitive assay of 17-hydroxyprogesterone.**

以下为常规应用缩写的中文注解

WB: 免疫印迹

ICC/IF: 细胞染色

IHC-P: 石蜡切片

IHC-Fr: 冰冻切片

## 属性

形式	Liquid
存储溶液	PBS
保存剂	0.09% Sodium azide
存放说明	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C.
浓度	Batch dependent (Please refer to the vial label for the specific concentration.)
偶联	Unconjugated



For full product information, images and publications, please visit our [website](#).

## 注意事项

仅供实验室使用。不适用于人类或动物的任何临床、治疗或诊断用途。不适合动物或人类食用。



For full product information, images and publications, please visit our [website](#).