

NCC antibody

Cat. No. GTX41969

宿主	Rabbit
克隆	Polyclonal
同种型	IgG
实验应用	WB, ICC/IF, IHC-P
种属反应	Human, Mouse, Rat

引用文献 (1)

实验应用


应用说明

We recommend a working dilution of 1:1000 (WB).

属性

形式	Liquid
存储溶液	0.1% BSA, 50% Glycerol
保存剂	1mM Sodium azide
存放说明	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
浓度	Batch dependent (Please refer to the vial label for the specific concentration.)
偶联	Unconjugated
注意事项	仅供实验室使用。不适用于人类或动物的任何临床、治疗或诊断用途。不适合动物或人类食用。

产品图片




GTX41969 ICC/IF Image

IF analysis of NCC in rat kidney tissues, showing using a 1:200 dilution of GTX41969. DAPI merge.



For full product information, images and publications, please visit our [website](#).

Date 2025 / 06 / 30 Page 1 of 2



GTX41969 WB Image

Western blot analysis of NCC in rat tissues using a 1:1000 dilution of GTX41969.



For full product information, images and publications, please visit our [website](#).