

Human Interferon omega protein (active)

Cat. No. GTX48194-pro

实验应用 Functional Assay

Speciess Human

Package

100 µg

实验应用

应用说明

The ED₅₀ was determined by a cytotoxicity assay using Human TF-1 cells is ≤ 0.01 ng/ml, corresponding to a specific activity of $\geq 1 \times 10^6$ units/mg.

属性

形式 Lyophilized powder

存储溶液 Batch dependent (please contact us for details)

存放说明 Store at -20°C or below. After reconstitution, keep as concentrated solution. Aliquot and avoid freeze-thaw cycles.

Region/Sequence MCDLPQNHGL LSRNTLVLLH QMRRISPFLC LKDRRDFRFP QEMVKGSQLQ KAHVMSVLHE MLQQIFSLFH TERSAAWNM
TLLDQLHTGL HQQLQHLETC LLQVVGEGES AGAISSPALT LRRYFQGIRV YLKEKKYSDC AWEVVRMEIM KSLFLSTNMQ
ERLRSKDRDL GSS

表达系统 E. coli

纯度 > 98% by SDS-PAGE and HPLC.

内毒素 < 1 EU/µg

偶联 Unconjugated

注意事项 仅供实验室使用。不适用于人类或动物的任何临床、治疗或诊断用途。不适合动物或人类食用。



For full product information, images and publications, please visit our [website](#).