

## Her6 blocking peptide

Cat. No. GTX48518-PEP

**Species** Zebrafish**Package**  
100 µg

## 实验应用

**产品说明** The peptide was used in the production of GTX48518.

## 属性

**形式** Lyophilized powder**存储溶液** Reconstitute with 200µl distilled water to obtain a 0.5mg/ml peptide solution. Lyophilized from 5% (v/v) acetonitrile/H<sub>2</sub>O.**保存剂** No preservatives**存放说明** Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.**Region/Sequence** KRKIKQTVTSLPHQ-C**注意事项** 仅供实验室使用。不适用于人类或动物的任何临床、治疗或诊断用途。不适合动物或人类食用。For full product information, images and publications, please visit our [website](#).