

SARS-CoV-2 Neutralizing Antibody ELISA Kit (Omicron BA.4 / BA.5 / BF.7 / BQ.1 / XBB.1.5)

Cat. No. GTX538288

实验应用 Neutralizing/Inhibition

种属反应 SARS Coronavirus 2

引用文献 (2)

Package

96 test

产品说明

摘要

The GeneTex SARS-CoV-2 Neutralizing Antibody ELISA Kit (Omicron BA.4/5 / BF.7 / BQ.1 / XBB1.5) is an in vitro assay for qualitative SARS-CoV-2 NAb screening. The kit includes recombinant His-tagged SARS-CoV-2 Spike RBD proteins for the Omicron variant BA.4/5, Omicron variant BF.7, Omicron variant BQ.1, and Omicron variant XBB 1.5, for Omicron variant-specific neutralization analysis.


属性

存放说明 The entire kit can be stored at 2-8°C for up to 6 months.

偶联 Horseradish peroxidase(HRP)


注意事项 仅供实验室使用。不适用于人类或动物的任何临床、治疗或诊断用途。不适合动物或人类食用。

产品图片



GTX538288 Neutralizing/Inhibition Image

Inhibition analysis of SARS-CoV-2 (COVID-19) Spike RBD antibody [HL1004] (GTX635793) (104-13333 pM) using SARS-CoV-2 Neutralizing Antibody ELISA Kit (Omicron BA.4 / BA.5 / BF.7 / BQ.1 / XBB.1.5).




GTX538288 Neutralizing/Inhibition Image

Inhibition analysis of SARS-CoV-2 (COVID-19) Spike RBD Omicron antibody [HL1866] (GTX637591) (104-13333 pM) using SARS-CoV-2 Neutralizing Antibody ELISA Kit (Omicron BA.4 / BA.5 / BF.7 / BQ.1 / XBB.1.5).



For full product information, images and publications, please visit our [website](https://www.genetex.com).



GTX538288 Neutralizing/Inhibition Image



For full product information, images and publications, please visit our [website](#).