

Human CCDC23 protein, His tag

Cat. No. GTX57419-pro

Species Human

Package
100 µg

实验应用

Observed MW 8 kDa.

属性

形式 Liquid

存储溶液 PBS, 20% Glycerol, 1mM DTT

保存剂 No preservative

存放说明 Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

浓度 0.5 mg/ml (Please refer to the vial label for the specific concentration.)

Region/Sequence Full length protein, N-terminal His-Tag; MGSSHHHHHH SSGLVPRGSH MGSEFMDPPA RKEKTKVKES VSRVEKAKQK
SAQQLKQRQ RAEIYALNRV MTELEQQQFD EFCKQMPPG E


表达系统 E. coli

纯度 > 90% by SDS-PAGE.

偶联 Unconjugated

注意事项 仅供实验室使用。不适用于人类或动物的任何临床、治疗或诊断用途。不适合动物或人类食用。

产品图片



GTX57419-pro Image

3µg Human CCDC23 protein (GTX57419-pro) by SDS-PAGE under reducing condition and visualized by coomassie blue stain.



For full product information, images and publications, please visit our [website](https://www.genetex.com).