

E. coli groL protein

Cat. No. GTX57497-pro

Speciess E. coli

Package
100 µg

实验应用

Observed MW 57 kDa.

属性

形式 Liquid

存储溶液 25mM Tris-HCl, 100mM NaCl, 10% Glycerol, 5mM DTT

保存剂 No preservatives

存放说明 Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

浓度 1 mg/ml (Please refer to the vial label for the specific concentration.)

Region/Sequence

Full length protein, MAAKDVKFGN DARVKMLRGV NVLADAVKVT LGPKGRNVWL DKSFGAPTIT KDGVSVAREI ELEDKFENMG
AQMVKVASK ANDAAGDGT TATVLAQAI TEGLKAVAAG MNPMDLKRG DKAVTAAVEE LKALSVPCSD SKAIAQVGTI
SANSDETGVK LIAEAMDKVG KEGVITVEDG TGLQDELDVV EGMQFDRGYL SPYFINKPET GAVELESPFI LLADKKISNI REMLPVLEAV
AKAGKPLII AEDVEGEALA TLVNTMRGI VKVAAPKAPG FGDRRKAMLQ DIATLTGGTV ISEEIGMELE KATLEDLGQA
KRVINKDIT TIIDGVGEEA AIQGRVAQIR QQIEEATSDY DREKLQERVA KLAGGVAVIK VGAATEVEMK EKARVEDAL
HATRAAVEEG VVAGGGVALI RVASKLADLR GQNEQDNVGI KVALRAMEAP LRQIVLNCGE EPSVVANTVK GGDGNYGYNA
ATEEYGNMID MGILDPTKVT RSALQYAASV AGLMITTECM VTDLPKNDA DLGAAGGMGG MGGMGGMM

表达系统 E. coli


纯度 > 95% by SDS-PAGE.

偶联 Unconjugated

注意事项 仅供实验室使用。不适用于人类或动物的任何临床、治疗或诊断用途。不适合动物或人类食用。

For full product information, images and publications, please visit our [website](#).

產品圖片



GTX57497-pro Image

3μg groL protein (GTX57497-pro) by SDS-PAGE under reducing condition and visualized by coomassie blue stain.



For full product information, images and publications, please visit our [website](#).