

GLP1 antibody [10]

Cat. No. GTX60952

宿主	Mouse
克隆	Monoclonal
同种型	IgG2a
实验应用	WB, ELISA
种属反应	Human

实验应用

应用说明

*最佳稀释倍数与浓度应由研究人员确认

Suggested dilution	Recommended dilution
WB	Assay dependent
ELISA	Assay dependent

以下为常规应用缩写的中文注解

WB: 免疫印迹

ICC/IF: 细胞染色

IHC-P: 石蜡切片

IHC-Fr: 冰冻切片


属性

形式	Liquid
存储溶液	10mM Phosphate, 500mM NaCl
保存剂	15mM Sodium azide
存放说明	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C.
浓度	Batch dependent (Please refer to the vial label for the specific concentration.)
偶联	Unconjugated
注意事项	仅供实验室使用。不适用于人类或动物的任何临床, 治疗或诊断用途。不适合动物或人类食用。



For full product information, images and publications, please visit our [website](#).

產品圖片

**GTX60952 ELISA Image**

It is possible to measure intact exendin-4 by a sandwich ELISA which is not subject to interference from GLP-1 or other glucagon-related peptides that may be present. GTX60952 is used as capture antibody in combination with a biotinylated detection antibody GTX60993-02.



For full product information, images and publications, please visit our [website](#).