

Mouse GI24 protein, human IgG1 Fc tag

Cat. No. GTX65677-pro**Species** Mouse**Package**

50 µg

属性

形式 Lyophilized powder**存储溶液** Reconstitute with 50µl ddH₂O to 1mg/ml. Lyophilized from 0.2µm-filtered PBS.**保存剂** No preservatives**存放说明** Store at -20°C or below. After reconstitution, keep as concentrated solution. Aliquot and avoid freeze-thaw cycles.**Region/Sequence** The extracellular domain of mouse VISTA (a.a. 33-191) is fused at the C-terminus to the Fc portion of human IgG1.**表达系统** HEK293 cells**纯度** ≥95% by SDS-PAGE.**内毒素** < 0.01 EU/µg**偶联** Unconjugated**注意事项** For In vitro laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.For full product information, images and publications, please visit our [website](#).