

PMCA antibody [5F10]

Cat. No. GTX72399

宿主	Mouse
克隆	Monoclonal
同种型	IgG1
实验应用	WB
种属反应	Human

引用文献 (1)

实验应用

应用说明

For WB: Use at a dilution of 1:200. Optimal dilutions/concentrations should be determined by the researcher.

属性

形式	Liquid
存储溶液	PBS, 0.2% BSA
保存剂	0.1% Sodium azide
存放说明	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. DO NOT FREEZE.
浓度	0.2 mg/ml (Please refer to the vial label for the specific concentration.)
偶联	Unconjugated
RRID	AB_383600
注意事项	仅供实验室使用。不适用于人类或动物的任何临床，治疗或诊断用途。不适合动物或人类食用。



For full product information, images and publications, please visit our [website](#).