

## Goat Anti-Rat IgG antibody, pre-adsorbed (HRP)

Cat. No. GTX73151

宿主	Goat
克隆	Polyclonal
同种型	IgG
实验应用	WB, IHC-P, IHC-Fr, ELISA
种属反应	Rat

引用文献 (3)  
Package  
500 µg

## 实验应用

## 应用说明

\*最佳稀释倍数与浓度应由研究人员确认

Suggested dilution	Recommended dilution
WB	Assay dependent
IHC-P	Assay dependent
IHC-Fr	1/50
ELISA	1/1000

以下为常规应用缩写的中文注解

WB: 免疫印迹

ICC/IF: 细胞染色

IHC-P: 石蜡切片

IHC-Fr: 冰冻切片

产品说明 Cross-reactivity with mouse immunoglobulins has been removed by solid phase adsorption.

## 属性

形式	Liquid
存储溶液	PBS, HRP stabilizer
保存剂	0.01% Thimerosal
存放说明	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. DO NOT FREEZE.
浓度	0.5 mg/ml (Please refer to the vial label for the specific concentration.)
免疫原	Rat IgG
纯化方式	Purified by affinity chromatography
偶联	Horseradish peroxidase(HRP)
RRID	AB_809601



For full product information, images and publications, please visit our [website](#).

**注意事项**

仅供实验室使用。不适用于人类或动物的任何临床, 治疗或诊断用途。不适合动物或人类食用。



For full product information, images and publications, please visit our [website](#).