

## CD2 antibody [MRC OX-34]

Cat. No. GTX75123

宿主	Mouse
克隆	Monoclonal
同种型	IgG2a
实验应用	IHC-P, IHC-Fr, FCM, IP
种属反应	Rat

## 实验应用

## 应用说明

For FACS: Use at a dilution of 1:25 - 1:75. For IHC-Fr, IHC-P and IP: Use at an assay dependent dilution. Optimal concentrations/dilutions should be determined by the researcher.

## 属性

形式	Liquid
存储溶液	PBS
保存剂	0.02% Sodium azide
存放说明	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
浓度	1.0 mg/ml (Please refer to the vial label for the specific concentration.)
偶联	Unconjugated
RRID	AB_379119
注意事项	仅供实验室使用。不适用于人类或动物的任何临床, 治疗或诊断用途。不适合动物或人类食用。



For full product information, images and publications, please visit our [website](#).