

## Integrin alpha 2 antibody [AK7]

Cat. No. GTX76461

宿主	Mouse
克隆	Monoclonal
同种型	IgG1
实验应用	IHC-Fr, FCM, Neutralizing/Inhibition
种属反应	Human, Baboon, Cynomolgus monkey, Rhesus Monkey

## 实验应用

## 应用说明

\*最佳稀释倍数与浓度应由研究人员确认

Suggested dilution	Recommended dilution
IHC-Fr	1/500-1/1,000
FCM	1/50-1/100
Neutralizing/Inhibition	Assay dependent

**Note : The epitope recognised by this antibody is reported to be sensitive to formaldehyde fixation and tissue processing. Recommend the use of acetone fixation for frozen sections.**

**Use 10µl of the suggested working dilution to label 10<sup>6</sup> cells in 100µl.**

**Mouse anti Human CD49b antibody, clone AK7 inhibits cell attachment to collagen.**

以下为常规应用缩写的中文注解

WB: 免疫印迹

ICC/IF: 细胞染色

IHC-P: 石蜡切片

IHC-Fr: 冰冻切片

## 属性

形式	Liquid
存储溶液	PBS
保存剂	0.09% Sodium azide
存放说明	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
浓度	1.0 mg/ml (Please refer to the vial label for the specific concentration.)
偶联	Unconjugated
RRID	AB_377368
注意事项	仅供实验室使用。不适用于人类或动物的任何临床、治疗或诊断用途。不适合动物或人类食用。



For full product information, images and publications, please visit our [website](#).