

Plasminogen antibody [4G6-B11]

Cat. No. GTX79690

宿主	Mouse
克隆	Monoclonal
同种型	IgG1
实验应用	WB, ELISA
种属反应	Human

实验应用

应用说明

*最佳稀释倍数与浓度应由研究人员确认

Suggested dilution	Recommended dilution
WB	Assay dependent
ELISA	Assay dependent

以下为常规应用缩写的中文注解

WB: 免疫印迹

ICC/IF: 细胞染色

IHC-P: 石蜡切片

IHC-Fr: 冰冻切片

属性

形式	Liquid
存储溶液	50mM Sodium Phosphate, 100mM NaCl, 1mM EDTA
保存剂	No preservative
存放说明	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
浓度	10 mg/ml (Please refer to the vial label for the specific concentration.)
偶联	Unconjugated
RRID	AB_11172358
注意事项	仅供实验室使用。不适用于人类或动物的任何临床、治疗或诊断用途。不适合动物或人类食用。



For full product information, images and publications, please visit our [website](#).

產品圖片



GTX79690 WB Image

WB analysis of various sample lysates using GTX79690 Plasminogen antibody [4G6-B11].

Lane 1 : Plasminogen Non-Reduced

Lane 2 : Plasmin Non-Reduced

Lane 3 : Kringle-I Fusion

Lane 4 : Plasminogen Reduced

Lane 5 : Plasmin Reduced

Lane 6 : M\W Markers



For full product information, images and publications, please visit our [website](#).