

## PRR7 antibody [TRAP3/10]

Cat. No. GTX79873

宿主	Mouse
克隆	Monoclonal
同种型	IgG2a
实验应用	WB, ICC/IF
种属反应	Human, Mouse, Rat

## 实验应用

## 应用说明

\*最佳稀释倍数与浓度应由研究人员确认

Suggested dilution	Recommended dilution
WB	Assay dependent
ICC/IF	Assay dependent

以下为常规应用缩写的中文注解

WB: 免疫印迹

ICC/IF: 细胞染色

IHC-P: 石蜡切片

IHC-Fr: 冰冻切片


## 属性

形式	Liquid
存储溶液	PBS
保存剂	15mM Sodium azide
存放说明	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. DO NOT FREEZE.
浓度	1 mg/ml (Please refer to the vial label for the specific concentration.)
偶联	Unconjugated
RRID	AB_11179259
注意事项	仅供实验室使用。不适用于人类或动物的任何临床, 治疗或诊断用途。不适合动物或人类食用。




For full product information, images and publications, please visit our [website](#).

## 產品圖片

**GTX79873 WB Image**

WB analysis of mouse brain tissue lysate using GTX79873 PRR7 antibody [TRAP3/10].

**GTX79873 ICC/IF Image**

ICC/IF analysis of HEK-293 cells co-transfected with PRR7 and GFP-PSD95 using GTX79873 PRR7 antibody [TRAP3/10].

Green : PSD95

Red : Primary antibody

Blue : DAPI



For full product information, images and publications, please visit our [website](#).