

Cytokeratin 3 antibody [AE5]

Cat. No. GTX80319

宿主	Mouse
克隆	Monoclonal
同种型	IgG
实验应用	WB, IHC-P
种属反应	Human, Rabbit, Bovine

实验应用

应用说明

Recommended Starting Dilutions:

For WB: Use at 1:1000-1:3000

Optimal dilutions should be determined experimentally by the researcher.

属性

形式	Liquid
存储溶液	PBS
保存剂	0.1% Sodium azide
存放说明	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
浓度	Batch dependent (Please refer to the vial label for the specific concentration.)
偶联	Unconjugated
注意事项	仅供实验室使用。不适用于人类或动物的任何临床, 治疗或诊断用途。不适合动物或人类食用。



For full product information, images and publications, please visit our [website](#).