

## CGRP antibody [4901]

Cat. No. GTX82726

宿主	Mouse
克隆	Monoclonal
同种型	IgG2a
实验应用	WB, ICC/IF, IHC-Fr, ELISA, IHC, Neutralizing /Inhibition, RIA
种属反应	Human, Mouse, Rat, Dog

引用文献 ( 10 )

## 实验应用


## 应用说明

Optimal working dilutions should be determined experimentally by the end user.

## 属性

形式	Liquid
存储溶液	10mM HEPES, 150mM NaCl, 0.01% BSA, 50% Glycerol
保存剂	No preservatives
存放说明	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
浓度	Batch dependent (Please refer to the vial label for the specific concentration.)
偶联	Unconjugated
RRID	AB_11173100
注意事项	仅供实验室使用。不适用于人类或动物的任何临床，治疗或诊断用途。不适合动物或人类食用。

## 产品图片



## GTX82726 ICC/IF Image

Immunostaining of CGRP (red) in cryopreserved DRG neurons grown on coverslips. 40xmag. Green staining is PGP9.5. Yellow staining is the colocalization of CGRP and PGP9.5.



For full product information, images and publications, please visit our [website](#).