

Cytokeratin 19 antibody [CK 202]

Cat. No. GTX82955

宿主	Mouse
克隆	Monoclonal
同种型	IgG1
实验应用	ICC/IF, IHC
种属反应	Human

实验应用

应用说明

IHC, ICC (frozen or formalin-fixed paraffin-embedded (FFPE) tissue sections, cell smears) For IHC dilute conc. antibodies 1:50-1:100, use streptavidin~biotin system or polymer system, incubate 30 minutes at room temperature.

属性

形式	Liquid
存储溶液	PBS
保存剂	0.05% Sodium azide
存放说明	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. DO NOT FREEZE.
浓度	0.20 mg/ml (Please refer to the vial label for the specific concentration.)
偶联	Unconjugated
注意事项	仅供实验室使用。不适用于人类或动物的任何临床, 治疗或诊断用途。不适合动物或人类食用。



For full product information, images and publications, please visit our [website](#).