

## NOVA1 blocking peptide

Cat. No. GTX88815-PEP

Species	Human
---------	-------

Package	100 µg
---------	--------

## 实验应用

产品说明	The peptide was used in the production of GTX88815.
------	---

## 属性

形式	Lyophilized powder
----	--------------------

存储溶液	Reconstitute with 200µl distilled water to obtain a 0.5mg/ml peptide solution. Lyophilized from 5% (v/v) acetonitrile/H <sub>2</sub> O.
------	---

保存剂	No preservatives
-----	------------------

存放说明	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
------	--

Region/Sequence	C-REMPQNVAKTEPVS
-----------------	------------------

注意事项	仅供实验室使用。不适用于人类或动物的任何临床、治疗或诊断用途。不适合动物或人类食用。
------	--



For full product information, images and publications, please visit our [website](#).