

Hamartin / TSC1 (isoform 1) blocking peptide

Cat. No. GTX89366-PEP

Species Human**Package**
100 µg

实验应用

产品说明 The peptide was used in the production of GTX89366.

属性

形式 Lyophilized powder**存储溶液** Reconstitute with 200µl distilled water to obtain a 0.5mg/ml peptide solution. Lyophilized from 5% (v/v) acetonitrile/H₂O.**保存剂** No preservatives**存放说明** Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.**Region/Sequence** C-QLHIMDYNETHHEH**注意事项** 仅供实验室使用。不适用于人类或动物的任何临床、治疗或诊断用途。不适合动物或人类食用。For full product information, images and publications, please visit our [website](#).